Risinno X1, chytrý prsten

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Risinno X1, Smart Ring

User Manual

Please read this manual carefully before using the device for the first time.
Please keep this manual for future reference.

The images and illustrations shown may vary depending on the device model.

1. Charge the ring.


Google Play

App Store

2. Find and install the **JYouPro** app.

3. Select your ring in the app and pair it.

4. A red flashing light means the ring is charging.

4. A green flashing light means the ring is fully charged

Case, built-in 200 mAh battery.

JYouPro

Connecting the ring

Open the **JYouPro** app, register an account, find the device, scan the device to locate the corresponding Bluetooth name, click to connect, and once connected, the ring's information will be displayed.

Personal Settings

In the personal information interface, you can configure related settings, step goals, units, and other information.

Sports

In the mobile app interface, go to the exercise section; you can select various exercise modes, including: walking, running, horseback riding, climbing, swimming, etc. Tap the exercise button to start tracking distance, pace, duration, and heart rate data for your activity. (This feature requires GPS to be enabled on your phone.

Activity Data

On the mobile app's home screen, you can view the number of steps, distance, duration, and calories burned recorded that day. You can also tap to view your fitness data for the past week or month. The number of step goals can be set in the app's personal information interface.

Heart Rate

On the mobile app's home screen, select "Heart Rate" heart rate, tap to start the measurement; the indicator light will turn on, and the measurement will last 30 to 60 seconds.

Blood Oxygen

On the mobile app's home screen, select "Blood Oxygen," tap the button to start the measurement; the indicator will light up green; continue testing for 30–60 seconds to complete the measurement, and the mobile app will display your current blood oxygen level, blood oxygen saturation, and fatigue level.

Sleep

The watch can display your total sleep time from the previous night and the duration of light and deep sleep. You can also tap to view sleep data for the past week or month.

Common Issues

1. Causes of Abnormal Heart Rate and Blood Oxygen Data Measured by the Ring

- 1) This is generally caused by poor contact between the ring's sensor and the body during measurement.
- 2) When taking a measurement, make sure the sensor is in full contact with your finger.

2. Causes of abnormal sleep data

- 1) Sleep monitoring is designed to simulate the state that people themselves perceive as the time they fall asleep and wake up; it must be worn correctly.
- 2) Going to sleep too late or sleeping before putting on the device may result in errors.
- 3) Sleep data is not tracked during the day, and the default sleep tracking period is from 9:30 PM to 12:00 PM the following day.

For more frequently asked questions, see the [My > FAQ] section in the app.

Announcement

1. The measurement results for this product are for reference only. Do not use them for medical purposes or as a basis for a doctor's instructions. The measurement results are intended for personal diagnosis.

2. This product has an IP68 waterproof rating and should not be used for deep-sea diving or prolonged immersion in water.

3. We reserve the right to change the content of this manual without prior notice. Some features may vary depending on the software version, which is normal.

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